



SPECIFICATIONS

VOLUME	300 ml
MAX PRESSURE	400 psi / 28 bar
INNER ID	54.00 mm
SCALE LENGTH	130.99 mm
PISTON STROKE	190 mm
\varnothing x LENGTH	69 x 211 mm
WEIGHT	2265 g / 4.99 lb




SPECIFICATIONS

VOLUME	200 ml
MAX PRESSURE	500 psi / 35 bar
INNER ID	44.00 mm
SCALE LENGTH	131.53 mm
PISTON STROKE	190 mm
Ø x LENGTH	59 x 211 mm
WEIGHT	1825 g / 4.02 lb




SPECIFICATIONS

VOLUME	100 ml
MAX PRESSURE	700 psi / 48 bar
INNER ID	34.90 mm
SCALE LENGTH	104.53 mm
PISTON STROKE	173 mm
Ø x LENGTH	50 x 193 mm
WEIGHT	1125 g / 2.48 lb




SPECIFICATIONS

VOLUME	50 ml
MAX PRESSURE	1000 psi / 69 bar
INNER ID	28.60 mm
SCALE LENGTH	77.83 mm
PISTON STROKE	151 mm
\varnothing x LENGTH	42 x 171 mm
WEIGHT	815 g / 1.8 lb




SPECIFICATIONS

VOLUME	20 ml
MAX PRESSURE	2300 psi / 159 bar
INNER ID	19.13 mm
SCALE LENGTH	69.58 mm
PISTON STROKE	151 mm
\varnothing x LENGTH	33 x 165 mm
WEIGHT	545 g / 1.2 lb



SPECIFICATIONS

VOLUME	10 ml
MAX PRESSURE	4000 psi / 276 bar
INNER ID	12.6 mm
SCALE LENGTH	80.02 mm
PISTON STROKE	151 mm
Ø x LENGTH	26 x 165 mm
WEIGHT	290 g / 0.64 lb



SPECIFICATIONS

VOLUME	6 ml
MAX PRESSURE	5000 psi / 345 bar
INNER ID	9.53 mm
SCALE LENGTH	84.11 mm
PISTON STROKE	151 mm
Ø x LENGTH	22 x 165 mm
WEIGHT	210 g / 0.46 lb